



Designed by the US military



(TACTICAL MECHANICAL TOURNIQUET)

Extremity hemorrhage is a major cause of preventable death globally and the number-one cause of preventable death on the battlefield.<sup>1</sup> Death from blood loss can occur in as little as five minutes from a severe extremity injury. Tourniquets are a proven device to stop severe, life-threatening bleeding.

Developed by the US military and recommended by the Committee on Tactical Combat Casualty Care (CoTCCC), the TMT™ is used by militaries, medical professionals, and public safety agencies worldwide.



**NO INTERNAL BAND:** Eases slack removal and allows for rapid resetting.

**2" WIDE:** Achieves occlusion at lower pressures.

**UV-PROTECTIVE ADDITIVE TO POLYMERS:** Durable enough to be stored outside of packaging.

**TWO LOCKING MECHANISMS:** Double-routed speed buckle and hook and loop backed strap prevents slippage and pressure loss.

**QUICK-DISCONNECT BUCKLE:** Prevents misrouting of tail and allows for trapped-limb application.

**ENHANCED TORSION BAR SECURING:** Audible signature "click" identifies the torsion bar is secure.

**MULTIFUNCTION:** Undeniable advantage when combined with the TQ-RAM† for point pressure over wound packing and additionally provides proximal hemorrhage control.

**SPECIFICATIONS**

PN	PDF185 (English IFU) PDF187 (Multi-language IFU)	31-115 (Berry Compliant)
NSN	n/a	6515-01-656-6191
CE marked	✓	✗
Purchase	Direct	DAPA, GSA, VA DAPA
Misc.	Accommodates limb circumferences 5.75" – 38" (14.6cm to 96.5cm)	Accommodates limb circumferences 5.75" – 38" (14.6cm to 96.5cm)
Packaged dimensions	4.5" x 2.5" x 2" (11.4 x 6.4 x 5.1cm)	4.5" x 2.5" x 2" (11.4 x 6.4 x 5.1cm)
Weight	3.5oz (99.2g)	3.5oz (99.2g)
Shelf life*	7 years	10 years

\* From date of manufacture

1. Eastridge, Brian J et al. "Death on the battlefield (2001-2011): Implications for the Future of Combat Casualty Care." *The Journal of Trauma and Acute Care Surgery* vol. 73,6 Suppl 5 (2012): S431-7. doi:10.1097/TA.0b013e3182755dcc

† May not be available in your region



visit [safeguardmedical.com](http://safeguardmedical.com)  
call (US): +1 (855) 428-6074 / (Int'l): +44 (0) 1568 613942  
email (US): [info@safeguardmedical.com](mailto:info@safeguardmedical.com) / (Int'l): [sales.admin@safeguardmedical.com](mailto:sales.admin@safeguardmedical.com)

